


The gold standard for measuring direct, real-time complete cardiac function.

Utilize the Pressure-Volume Analysis Module for LabChart to measure and analyze in real-time or post-acquisition *in vivo* ventricular pressure-volume in small and large animals, or *ex vivo* using working heart systems.

This Module offers smart presets for different animals and streamlined workflows guiding you step-by-step from calibration through to analysis. When combined with Millar's extensive range of PV catheters, you have a highly sensitive, minimally invasive way to measure and analyze changes in ventricular pressure-volume in small and large animals. PV Loops are simultaneously generated in Loop View as pressure and volume are being recorded.

## Applications

- Systolic and diastolic dysfunction
- Valvular heart disease
- Cardiomyopathy
- Pulmonary hypertension
- Cardiac failure
- Cardiovascular remodeling and vascular occlusions
- Ischemia / Reperfusion studies
- Cardiac resynchronization therapy



Above: Loop View (right) generated from pressure and volume data in Chart View (left) provides online and offline analysis of all or selected PV Loops with ESPVR and EDPVR values.

# Powerful Analysis Options

This module provides a variety of automatically calculated hemodynamic data sets with several analysis options.

## Streamlined Workflow

The Pressure-Volume Analysis Module guides your journey from pressure, conductance, and volume calibration through to data analysis with supporting workflows for large animals, small animals, or echocardiogram users.

## Loop View

Display and analyze individual or selected loops of interest. Calculate and display end-systolic and end-diastolic PV relationships. Apply linear, exponential, or quadratic regression fit options to your ESVPR and EDVPR data.

## Linear regression information

Display your linear regression information, whether you're online or offline, with our built-in plots:

- Preload-recrutable stroke work (PRSW)
- Preload-recrutable maximum pressure change (dP/dt Max versus EDV)
- Preload-recrutable total mechanical work (PVA versus EDV)
- Afterload-dependent total mechanical work (PVA versus ESP)

## Calibration of Data

Preset calibrations for small and large mammals enable you to follow best practices for recording and analyzing PV data using Millar catheters. The Workflow modes have been optimized for use with the MPVS Ultra Foundation Systems.

- Relative Ventricular Volume Calibration
- Cuvette Calibration for absolute volumes. This module includes known volumes of common cuvettes.
- Saline Calibration to correct for parallel conductance or parallel volume effects.
- Or simplify your calibration using data from an external echocardiogram

## Hemodynamics Table


Calculate and display 29 hemodynamic parameters (including pressure-volume area, stroke work, and cardiac output). Easily export your data to LabChart DataPad or as csv/text files.

*Below:* Hemodynamics Table

Pressure Volume Recording: Hemodynamics Table																												
SW (mmHg*RVU)	CO (RVU/min)	SV (RVU)	Vmax (RVU)	Vmin (RVU)	Ves (RVU)	Ved (RVU)	Pmax (mmHg)	Pmin (mmHg)	Pmean (mmHg)	Pdev (mmHg)	Pes (mmHg)	Ped (mmHg)	HR (bpm)	EF (%)	Ea (mmHg*RVU)	PVA (mmHg*RVU)	PE (mmHg*RVU)	CE	Tau (ms)									
1 214.0	1273	2.101	9.645	7.544	8.625	9.087	145.4	-3.028	48.50	148.4	144.8	1.436	606.1	22.11	68.90	372.5	158.5	0.5744	4.048									
2 194.4	1223	2.018	9.547	7.529	8.578	8.929	142.1	-2.872	46.56	144.9	141.9	1.186	606.1	21.58	70.32	348.7	154.3	0.5576	3.842									
3 202.9	1251	2.043	9.558	7.515	8.572	8.972	143.2	-3.059	47.87	146.3	142.7	1.030	612.2	21.77	69.82	358.6	155.7	0.5659	3.961									
4 213.5	1211	1.998	9.506	7.509	8.521	9.285	140.5	-2.580	46.20	143.0	140.2	1.218	606.1	21.06	70.18	364.0	150.5	0.5867	4.159									
5 204.0	1224	2.020	9.510	7.480	8.548	9.031	142.1	-3.081	46.67	145.2	141.6	1.249	606.1	21.45	70.13	355.7	151.7	0.5724	4.077									
6 176.5	1122	1.852	8.339	7.487	8.478	8.719	137.4	-2.810	44.16	140.2	137.3	0.5307	606.1	20.11	74.15	323.5	147.0	0.5457	4.112									
7 189.8	1170	1.931	9.416	7.485	8.484	8.843	140.9	-3.028	46.17	143.9	140.4	0.6244	606.1	20.89	72.68	341.9	152.1	0.5551	3.986									
8 192.5	1130	1.864	9.304	7.440	8.372	9.125	137.0	-2.488	44.38	139.5	136.0	1.499	606.1	20.07	72.94	338.5	146.0	0.5687	4.425									
9 181.7	1089	1.797	9.250	7.453	8.407	8.970	137.1	-3.028	43.92	140.2	136.9	1.611	606.1	19.38	76.16	325.6	143.9	0.5581	4.256									
10 157.5	990.8	1.635	9.093	7.458	8.297	8.680	133.1	-3.059	41.95	136.1	132.6	0.4058	606.1	18.13	81.14	301.3	143.8	0.5227	4.148									
11 171.9	1065	1.757	9.243	7.485	8.388	8.749	138.0	-2.903	45.02	140.9	137.3	0.9990	606.1	19.27	78.16	320.9	149.1	0.5355	4.200									
12 188.2	1114	1.819	9.243	7.424	8.372	9.114	135.9	-2.510	44.36	138.2	135.2	1.717	612.2	19.71	74.32	328.3	140.1	0.5703	4.403									
13 181.6	1067	1.760	9.192	7.432	8.379	9.057	135.6	-2.747	43.43	138.3	135.4	1.623	606.1	19.19	76.94	321.2	139.6	0.5654	4.270									
14 159.3	979.0	1.632	9.071	7.439	8.278	8.660	133.3	-2.903	42.00	136.2	132.9	0.4683	600.0	18.17	81.45	300.9	141.6	0.5293	4.167									
15 162.8	1002	1.653	9.099	7.447	8.313	8.698	135.6	-2.935	43.78	138.6	135.2	1.061	606.1	18.37	81.84	305.4	142.6	0.5330	4.059									
16 180.0	1073	1.752	9.138	7.386	8.275	9.027	132.7	-2.373	43.04	135.1	132.1	1.155	612.2	19.20	75.40	315.2	135.2	0.5712	4.392									
17 176.6	1028	1.713	9.120	7.407	8.292	9.093	131.4	-2.716	41.21	134.1	131.1	1.592	600.0	18.78	76.36	310.5	133.9	0.5687	4.478									
18 155.4	955.3	1.576	8.989	7.412	8.216	8.649	130.9	-3.122	41.40	134.0	130.6	0.9366	606.1	17.64	82.83	290.5	135.0	0.5351	4.211									
Avg	155.8	1.588	9.086	7.398	8.355	8.740	127.2	-2.824	40.29	130.1	126.8	0.9691	604.3	17.72	81.06	287.2	131.4	0.5389	4.293									
Min	95.37	634.8	1.058	8.379	7.270	7.738	8.126	-3.497	32.73	108.0	105.3	0.0000	600.0	12.69	68.90	211.6	106.7	0.4507	3.842									
Max	214.0	1273	2.101	9.645	7.544	8.625	145.4	-1.280	48.50	148.4	144.8	2.404	612.2	22.11	105.0	372.5	158.5	0.5867	4.760									
Count	47	47	47	47	47	47	47	47	47	47	47	47	47	47	47	47	47	47	47									

## LabChart Analysis Manager

LabChart allows you to analyze and compare separate data regions within a data file. Save the analysis results and settings within a single file. Quickly and easily reapply those saved analysis settings to other files.



*Left:* Analysis Manager

## Ordering Information

The PV Loop Module for LabChart can be purchased individually as an Add-On for LabChart 8 (MLS240/8 Windows only), or as part of LabChart Pro (MLS260/8).

LabChart Pro includes LabChart software and all LabChart Modules, providing powerful data acquisition and analysis capabilities.

PowerLab and LabChart are trademarks of ADInstruments Pty Ltd. All other trademarks are the property of their respective owners. Products supplied by ADInstruments are intended for use in research and teaching applications and environments only.



Visit [adinstruments.com](https://www.adinstruments.com) or contact your local ADInstruments representative for more information

Australia | Brazil | Europe | India | Japan | China | Middle East | New Zealand | North America | Pakistan | South America | South East Asia | United Kingdom

[adinstruments.com](https://www.adinstruments.com)

