

MLA1095 Electrode Paste (3 pack)

Accessories

MLA1095 Electrode Paste (3 pack) is a conductive, water-soluble adhesive paste for use with reusable EEG electrodes. When used with MLAWBT9 EEG Flat Electrodes, apply just enough paste inside the concave side of the electrode to slightly overfill it. Place the electrode on the prepared electrode site and press down with medium pressure.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1095 Electrode Paste (3 pack)

For use with:
MLAWBT9 EEG Flat Electrodes

Also Available: MLA1095B Electrode Paste (6 pack)