

## **MLA1210 Spring Clip Electrodes**

### *Accessories*

Set of 3 gold plated electrodes. Each electrode has a one meter lead wire ending in a 2mm pin. The electrodes are suitable for use with the FE136 Animal Bio Amp. They can also be used with the MLA2540 or MLA2340 Bio Amp Cables by using the MLAC14 Adapter Plugs. The electrodes are typically used for measuring electrical activity from the isolated heart in a Langendorff or Working Heart preparation.



#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA1210 Spring Clip Electrodes

For use with:

FE136 Animal Bio Amp

MLA2340 3 Lead Shielded Bio Amp Cable (with MLAC14 adapter pack)

MLA2540 5 Lead Shielded Bio Amp Cable (with MLAC14 adapter pack)