

MLA1214 Spring Clip Electrodes

Accessories

A set of three leads with 1.5 mm shrouded sockets at one end and gold plated spring clip electrodes at the other. It is designed for use with the FE136 Animal Bio Amp. The lead wires are 1.5 m in length. They are shipped non-sterile.

The electrodes are typically used for measuring electrical activity from the isolated heart in a Langendorff or Working Heart preparation.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1214 Spring Clip Electrodes

For use with:
FE136 Animal Bio Amp