

## MLA1505 Lead Wires (5 Alligator Clips)

### Accessories

A set of five alligator clip lead wires for use with ADInstruments Single and Dual Bio-Amps (via a Bio-Amp Cable), directly to an ADInstruments Quad or Octal Bio-Amp or directly to a GT201 16 Channel Bio-Amp. Suitable for use as a recording cable in experiments with the MLT012 Nerve Chamber.



### Caution

Read "Statement of Intended Use" on our website before use. Leads are not supplied sterile. Appropriate Infection Control Procedures should be followed between subsequent subject connections.

### Ordering Information:

MLA1505 Lead Wires (5 Alligator Clips)

For use with:


MLA1340 3 Lead Bio Amp Cable

MLA1540 5 Lead Bio Amp Cable

FE234 Quad Bio Amp

FE238 Octal Bio Amp

GT201 16 Channel Bio Amp



NOTE: If you have a Shielded Bio Amp Cable (MLA2340 & MLA2540). Lead Wires can be connected to the top row signal pins.