

MLA1515 Bio Cable (Dual Bio to 5 Alligator, 1m)

Accessories

The MLA1515 is a Bio Cable specifically designed for use in measuring biopotentials such as ECG and EMG from animals.

The cable has a 6-pin dual Bio connector to five Alligator clips and can be used with any Dual Bio Amp or PowerLab unit with a built in Dual Bio Amp.



Cable Colour	Description
Brown	Ch 1 – Negative
Green	Ch 1 – Positive
Black	Ch 2 – Negative
Red	Ch 2 – Positive
White	Reference/Common

Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1515 Bio Cable

For use with:

Any PowerLab with a built-in Dual Bio Amp

Any ADInstruments Dual Bio Amp