

MLA1610 Shielded Lead Wires (5 Micro-Hooks)

Accessories

A set of five Micro-Hook lead wires for use with MLA2540, MLA2350 or MLA2340 Shielded Bio Amp Cables. Suitable for use as a recording cable in experiments with the MLT016 Nerve Chamber.



Caution

Read “Statement of Intended Use” on our website before use. Leads are not supplied sterile. Appropriate Infection Control Procedures should be followed between subsequent subject connections.

Ordering Information:

MLA1610 Shielded Lead Wires (5 Micro-Hooks)

For use with:

MLA2340 3 Lead Shielded Bio Amp Cable

MLA2350 3 Lead Shielded Bio Amp Cable

MLA2540 5 Lead Shielded Bio Amp Cable

MLT016 Nerve Chamber