

## **MLA1615 Shielded Lead Wires (long Alligator Clip)**

### *Accessories*

A set of five Alligator Clip lead wires (1 m) for use with MLA2540, MLA2350 or MLA2340 Shielded Bio Amp Cables.



#### **Caution**

Read “Statement of Intended Use” on our website before use. Leads are not supplied sterile. Appropriate Infection Control Procedures should be followed between subsequent subject connections.

#### **Ordering Information:**

MLA1615 Shielded Lead Wires (long Alligator Clip)

For use with:

MLA2340 3 Lead Shielded Bio Amp Cable

MLA2350 3 Lead Shielded Bio Amp Cable

MLA2540 5 Lead Shielded Bio Amp Cable