

MLA2350 Shielded Bio Amp Cable (3 Lead)

Spare Parts

The MLA2350 Shielded Bio Amp Cable (3 Lead) is suitable for use with MLA2504 (3 snap-on) Shielded Lead Wires, MLA2505 Shielded Lead Wires (5 snap-on) and MLA1610 Shielded Lead Wires (5 Micro-Hooks).

The MLA2350 Shielded Bio Amp Cable (3 Lead) is suitable for use with the FE231 Bio Amp.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA2350 Shielded Bio Amp Cable (3-Lead)

For use with:

MLA2504 Shielded Lead Wires (3 snap-on)

MLA2505 Shielded Lead Wires (5 snap-on)

MLA1610 Shielded Lead Wires (5 Micro-Hooks)

FE231 Bio Amp