

MLA2505 Shielded Lead Wires (5 Snap-on)

Spare Parts

Shielded Lead Wires (MLA2505, 5 pack) are suitable for use with Shielded Bio Amp Cables (MLA2340, MLA2350 & MLA2540). The lead wires are 98 cm in length with 4 mm 'snap-on' connectors for use with Disposable ECG electrodes (MLA1010 & MLA1010B), MLA700 Reusable ECG Electrodes and MLA710 Chest ECG Electrodes. May be cold sterilized.



Caution

Read "Statement of Intended Use" on our website before use. Leads are not supplied sterile. Appropriate Infection Control Procedures should be followed between subsequent subject connections.

Ordering Information:

MLA2505 Shielded Lead Wires (5 snap-on)

For use with:

MLA2340 3 Lead Shielded Bio Amp Cable

MLA2540 5 Lead Shielded Bio Amp Cable

MLA1010 & MLA1010B Disposable ECG Electrodes

MLA2350 3 Lead Shielded Bio Amp Cable

MLA700 Reusable ECG Electrodes

MLA710 Chest ECG Electrodes