

## **MLA265 Stimulator Rod With Cable**

### *Accessories*

The Stimulating rod is 150 mm long with an OD of 10 mm. It has a rounded stainless steel tip that can be placed on the skin or tissue surface. The other end terminates in a stud connector. The rod is supplied with a cable 1800 mm long that terminates in 4 mm shrouded connectors and a press stud connector to connect to the rod.



It is suitable for use with PowerLabs with a built-in Stimulus Isolator or FE180 Stimulus Isolator.

#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### **Ordering Information:**

MLA265 Stimulator Rod with Cable

For use with:

FE180 Isolated Stimulator

Any PowerLab with built-in Stimulus Isolator