

MLA304 Disposable Filters

Accessories

A Pack of 50 Disposable Filters protect against cross contamination without impairing airflow. MLA1026 Reusable Mouthpieces and MLA1011A Clean Bore Tubing & Flow Head Adapter can be used with the filters. Intended for single human use only — do not wash or sterilize.



Dimensions:

Subject end 31 mm (1.2") outer diameter

Tubing end 34 mm (1.3") outer diameter

Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA304 Disposable Filters

For use with:

MLA1026 Reusable Mouthpieces

MLA1011A Clean Bore Tubing and Flow Head Adapter