

## Rodent Restrainers

### Accessories

Rodent restrainers supplied by ADInstruments are constructed of perspex. They are available in six sizes to suit different sized animals. Slots and holes are provided for ventilation and access for drug administration. Space inside the chamber can be adjusted using the slide piston.



#### Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

#### Ordering Information:

MLA5016 Rodent Restrainer for Mice up to 35 g

MLA5018 Rodent Restrainer for Mice up to 50 g

MLA5020 Rodent Restrainer for Rats (80 - 200 g)

MLA5022 Rodent Restrainer for Rats (180 - 320 g)

MLA5024 Rodent Restrainer for Rats (300 - 440 g)

