

MLA5030 Rodent Tail Cuff Holder

Accessories

The MLA5030 Rodent Tail Cuff Holder allows users of rodent restrainers to secure the rat or mouse tail transducer and cuff (MLT125M or MLT125R) when attached to a mouse or rat tail.

The Rodent Tail Cuff Holder is specifically manufactured for use with the ADInstruments supplied rodent restrainers: MLA5016, MLA5018, MLA5020, MLA5022, MLA5024, MLA5025.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA5030 Rodent Tail Cuff Holder

For use with:

MLA5016 Rodent Restrainer for Mice up to 35 g

MLA5018 Rodent Restrainer for Mice up to 50 g

MLA5020 Rodent Restrainer for Rats (80 - 200 g)

MLA5022 Rodent Restrainer for Rats (180 - 320 g)

MLA5024 Rodent Restrainer for Rats (300 - 440 g)

MLA5025 Rodent Restrainer for Rats (440 - 550 g)

