

# MLAIME EEG Clip Electrodes

## Accessories

The MLAIME EEG Clip Electrodes are a set of two ear clip earth electrodes, 9.5 mm (0.4") wide with 1.9 m (6.2') long black leads. The electrodes are fitted with shrouded 1.5 mm (0.06") sockets for connection to MLA2340 or MLA2540 Shielded Bio Amp Cables. These electrodes can also be used with the GT205 8 Channel, GT201 16 Channel Bio Amps and ML138 Octal Bio Amp.



### Caution


Read "Statement of Intended Use" on our website or in "Getting Started with PowerLab" before use. Clean with water and mild detergent. Do not use alcohol.

### Ordering Information:

MLAIME EEG Clip Electrodes

For use with:

- MLA2340 3 Lead Shielded Bio Amp Cable
- MLA2540 5 Lead Shielded Bio Amp Cable
- GT205 8 Channel Bio Amp
- GT201 16 Channel Bio Amp
- ML138 Octal Bio Amp



NOTE: Lead Wires are connected to the top row of signal pins in Shielded Bio Amp Cables