

## **SP3925 Membranes (10) and O-rings (8) for rat tail cuff**

### *Spare Parts*

A set of 10 replacement latex membranes for use on the pressure cuff of the rat MLT125/R pressure cuff and transducer of the rat NIBP System (FE125/R).



#### **Caution**

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

**Warranty is valid for 30 days from date of purchase and is void if the transducer is damaged due to neglect.**

#### **Ordering Information:**

SP3925 Membranes (10) and O-rings (8) for rat tail cuff

For use with:

MLT125/R pressure cuff