

MLA1028 Face Mask

Accessories

The MLA1028 Face Mask features a Y-shaped, two-way, non-rebreathing valve designed to separate inspired and expired gasses.

Typical applications include exercise physiology and respiratory teaching and research.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1028 Face Mask

For use with:

MLA1015 Breathing Tube

MLA245 Gas Mixing Chamber

MLA246 Gas Mixing Chamber