

MLA1029 Face Mask Kit (Adult, small & medium)

Accessories

The MLA1029 Face Mask Kit includes a head strap, Y-shaped two-way non-rebreathing valve, with interchangeable Adult Face Masks (medium or small). The unit is designed to separate inspired and expired air by using valves which allow flow in one direction only.

Typical applications include exercise physiology and respiratory teaching and research.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1029 Face Mask Kit Contains:

MLA1028 Face Mask Kit (Adult, medium)

SP0141 Adult Face Mask (small)

For use with:

MLA1015 Breathing Tube

MLA246 Gas Mixing Chamber