

MLA1081 Flow Head Adapter

Accessories

An adapter used to connect the MLA1011 Clean Bore Tubing to MLT300L, MLT1000L or MLT3813H Flow Heads.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1081 Flow Head Adapter

For use with:

MLT1000L Respiratory Flow Head 1000L

MLT300L Respiratory Flow Head 300L

MLT3813H Heated Pneumotach

MLA1011 Clean Bore Tubing