


MLA1204 Needle Electrodes (2 mm Pin)

Accessories

A pack of 5 needle electrodes with 2 mm (0.05") pin plugs for connection to the FE136 Animal Bio Amp.

Needle gauge:	29
Needle length:	12 mm (0.5")
Needle material:	Surgical steel
Lead wire:	122 cm (48")



Caution

Read "Statement of Intended Use" on our website or in "Getting Started with PowerLab" before use.

Ordering Information:

MLA1204 Needle Electrodes (2 mm Pin)

For use with:
FE136 Animal Bio Amp