

MLA1212 Micro-Hook Electrodes (3 pk)

Accessories

A set of three 1.5 mm shrouded sockets to Micro-Hook lead wires for use with the FE231 Bio Amp or older FE136 Animal Bio Amp. The lead wires are 75 cm in length with the Micro-Hooks constructed of molded nylon with a stainless steel spring and gold-plated beryllium copper hooks. A set is included with every FE231 Bio Amp.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1212 Micro-Hook Electrodes (3 pk)

For use with:

FE231 Bio Amp

FE136 Animal Bio Amp (obsolete)