


MLA1213 Needle Electrodes for FE231 (3 pk)

Accessories

A set of three 1.5 mm shrouded sockets to monopolar needle electrodes (29 gauge) lead wires for use with the FE231 Bio Amp or the obsolete FE136 Animal Bio Amp. The lead wires are 122 cm in length and the needle electrodes are made from surgical steel and 12 mm in length. They are shipped non-sterile.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA1213 Needle Electrodes for FE231 (3pk)

For use with:

FE231 Bio Amp

FE136 Animal Bio Amp