

MLA260 Stimulator Cable

Accessories

Cable with two shrouded 4 mm (0.16") connectors to two alligator clips suitable for use as a stimulator cable in experiments with the MLT012 Nerve Chamber. Connects to any Powerlab with built-in isolated Stimulator FE180 Stimulus Isolator.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA260 Stimulator Cable

For use with:

MLT012 Nerve Chamber

Any PowerLab with a built-in isolated stimulator

FE180 Stimulus Isolator