

MLA5530 Calibration Syringe 3 L

Accessories

The MLA5530 3 L Calibration Syringe provides an easy and reliable method of calibrating respiratory equipment. Other volumes ranging from 10 mL to 7 L are available on request.



When calibrating, a steady rate of movement of the plunger should be maintained. High rates of movement and sudden jolts may cause unwanted interference on the recorded signal.

Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLA5530 Calibration Syringe 3 L

For use with:

MLT1000L Flow Head