

MLA6024 Desiccant Cartridge

Accessories

A cartridge for removing moisture from a gas stream. Can be used multiple times and may be recharged by replacing the desiccant material with available refill (MLA6026).

Recommended for use with the ML206 Gas Analyser if analyzing gas from exercising subjects or in high humidity environments.

Size (L x B x H): 95 x 60 x 25 mm (3.7" x 2.3" x 0.98")



Caution

Read "Statement of Intended Use" on our website or in "Getting Started with PowerLab" before use.

Ordering Information:

MLA6024 Desiccant Cartridge

For use with:

MLA0343 Drying Tube

MLA6026 Desiccant Refill

ML206 Gas Analyzer