

## **MLA700 Reusable ECG Electrodes**

### *Accessories*

The MLA700 contains a set of three reusable clamp electrodes for use with MLA2503 or MLA2505 Shielded Lead Wires or MLA0313 and MLA0315 Lead Wires (Unshielded). The electrodes have a surface area of 24 cm<sup>2</sup> and can be attached to wrists or ankles. The electrodes provide good contact and low noise and are ideal for student use. May be cold sterilized only.



Recommended use with Electrode Cream (MLA1090) to minimize resistance in recordings.

#### **Caution**

Read “Statement of Intended Use” on our website before use. Leads are not supplied sterile. Appropriate Infection Control Procedures should be followed between subsequent subject connections.

#### **Ordering Information:**

MLA700 Reusable ECG Electrodes

For use with:

MLA2503 Shielded Lead Wires (3 snap on)

MLA2505 Shielded Lead Wires (5 snap on)

MLA1090 & MLA1090B Electrode Cream