

MLAC11 Grass Adapter Cable

Accessories

Adapter cable that connects a Grass FT03 or FT10 force transducer to the ML110 and ML221 Bridge Amp, ML112 and ML224 Quad Bridge Amp, ML228 Octal Bridge Amp.

NOTE: Not suitable for connecting to ML118 Quad Bridge Amp or ML119 Octal Bridge Amp. For these bridge amps use the MLAC18 Grass Adapter Cable.



Caution

Read “Statement of Intended Use” on our website or in “Getting Started with PowerLab” before use.

Ordering Information:

MLAC11 Grass Adapter Cable

For use with:

ML110 Bridge Amp

ML221 Bridge Amp

ML112 Quad Bridge Amp

ML224 Quad Bridge Amp

ML228 Octal Bridge Amp